

SKIN ANALYSIS SYSTEM

New professional system to analyze skin health and recommend cosmetic products and treatments



User friendly, specially designed to accomplish either professional and customers needs



Powerful customer management which offers professional and customized recommendations




Marketing tool to support sales and promotions of dermocosmetic products



Several applications that cater to the aesthetics centers, cosmetics shops, pharmacies, cosmetic agents, etc.



dermo prime SOFTWARE









TECHNICAL SPECIFICATIONS

Specially designed image capture technology, defines all skin issues which are then analyzed through the exclusive software DermoPrime. The results are based on mathematical algorithms, which are continuously improved according to a multimillion image database.





CAMERA

Wireless video transmission in real time, with high magnification images

 Image Sensor	5MEGA Pixel CMOS with Optic Zoom	
 Conection	Wifi, Wifi Direct	
 Weight	aprox. 150 g	
 Battery	1050 mAh (Li-Polymer), 60-90 min of Continuous Use · Fast charging conection	
 Image Mode	Skin	VSL Mode: Pores · Wrinkles XPL Mode: Spots · Keratin UVL Mode: Acne
		Integrated Hydration Sensor Sebum and Keratin indicators included

DOCKING / STANDER

Elegant and ergonomic design for organise and charge the camera

-  **Camera charger**
-  **Power supply connection**

dplite2



DIMENSIONS <i>(without the lens)</i>	29,5 (W) x 23 (D) x 36 (H) cm
MAX. CONSUMPTION	2,5 With standard lens
BATTERY	1050 mAh Li-Ion Polymer
WEIGHT <i>(without battery, with lens)</i>	130 g

MEASUREMENTS



MOISTURE



PORES



SPOTS



WRINKLES



ACNE



KERATIN



SKIN TYPE